

PHILIPS

Pagewriter TC10



A Portable Workflow

6 channel cardiograph – TC10

The Philips Pagewriter TC10 cardiograph is portable, affordable and easy to operate. It supports advanced Pagewriter cardiograph features, scalable network connectivity and World-class Philips DXL ECG algorithm.

The Intuitive 1-2-3 touch screen lets physician, technician, nurses and assistants efficiently acquire, analyze, store, print and access ECGs. With it's compact size and lightweight design, the Pagewriter TC10 is portable workflow solution that supports patient needs in a wide range of care settings

Key advantages

- Compact size
- Lightweight design
- Built in recorder
- Storage facility
- DXL algorithm for CA detection (option)
- ST analysis and ST MAP (option)
- Extended storage and PDF out (option)
- PDF for 12CH print
- Wireless networking option to iECG



Connect Leads

Check for loose leads.

Enter ID

Barcode ID to find record from worklist.

Take ECG

Print, save, transfer and retrieve previous ECG.



Features :

Signal Quality Indication (standard)

Lead Off advisory with lead map indication, lead colour to indicate the quality of individual leads, lead placement check, heart rate display, full screen preview, real time wave display with multiple choice of formats selection


DXL Algorithm and STEMI diagnostic aid (optional)

More than 600 interpretive statements, integrated pediatric analysis, morphology analysis in each of the 12 leads, ST analysis with MAP, CA detection and critical values

ECG Archive system (optional)

Archive with PDF / XML convertor for easy documentation, A4 size print facility with 12L view on a plain paper (PDF print via PC and printer)

Philips Healthcare is a part of Philips Electronics



Pagewriter TC10



*Clamps and Bulbs **



10 Lead patient cable

** Accessories shown in the brochure may no be a part of standard configuration*

