

## RESPIRATORY PRODUCTS

# AERO PISTON COMPRESSOR NEBULIZER



**Aero** nebulizer is designed for continuous use built with heavy duty oil free piston compressor. No Maintenance is required as the oil free piston compressor does not need any lubrication.

Through the **Aero** nebulizer, a treatment of 2cc of medicine is completed in 7 minutes, nebulized in very small particles measuring between 0.5 and 5.5 microns, reaching the deepest areas of the respiratory tract.

Dust and impurities are kept away by the filter which ensures a safe medical use & is easily replaceable.

Each **Aero** nebulizer is thoroughly tested and is covered with 1 year warranty.

# Aero

**PISTON COMPRESSOR NEBULIZER**

... Serving for Better Life Care