

## AERO OXYGEN CONCENTRATOR



### FEATURES

- + PSA (Pressure Swing Adsorption) Technology
- + Adjustable Flow rate
- + Large LCD Display for
  - Accumulating timing
  - Present timing
  - Presetting timing
- + Digital:
  - Pressure digital
  - Temperature digital
  - Purity digital
  - SPO2 digital
- + Maintenance Reminder

### SPECIFICATION

Flow rate	: 0-5 L/min
Purity :	: 93% ± 3% (MOLECULER SIEVE IMPORTED FROM USA)
Outlet Pressure (Mpa)	: 0.03 - 0.07
Sound Level (dB)	: < 50
Power	: AC220 V ±10%. 50/60 ±1 Hz
Power	: < 480 W
Timer	: Accumulating timing, Present timing
Alarm	: Power Failure, High & Low Pressure, High Temperature
Net Weight	: 30 kg.
Size (mm)	: 365 x 375 x 600