

## ELOXY 5W OXYGEN CONCENTRATOR



### FEATURES

- + Quiet Operation- less than 36Dba
- + Light weight. Easy to Transport and Move
- + Superior grade of molecular sieve
- + Exclusive Moisture reducing Wick Assembly
- + Spring mounted High Efficiency Thomas Nexus Compressor or equivalent
- + With inbuilt Nebulizer, Design so simple it does its job in best way

### SPECIFICATION

Electrical Requirement	: 230 Volt – 50/60Hz
Flow Delivery Rate-	: 0 to 5LPM
Oxygen Concentration-	: 93% + 3%
Power Consumption-	: 390 Watts
Operating Pressure (Mpa)-	: 0.04 – 0.07
Safety Alarms & Features-	: Power Failure
Filters-	: Cabinet, Compressor Intake & Bacteria
Battery-	: 9 Volt
Weight	: Model 5G – 18Kg, Model 5W – 16kg
Dimensions-	: 430mm L x 380mm W x 690mm H
Operating Environment	: Ambient Temperature: 50°F to 100°F (10°C to 40°C) Humidity: 15% to 95%, non – condensing
Storage Range	: Temperature: 0°F to 140°F (20°C to 60°C) Humidity: 15% to 95%, non – condensing