

## AERO OXYGEN CONCENTRATOR



### FEATURES

- + PSA - Pressure swing adsorption technology and adjustable flowrate.
- Nebulizer device can be fitted if necessary.
- Power failure alarm, pressure alarm, low purity alarm.
- Double flow meter control for two patients.
- Dehumidity function.
- Timer: (a) total working hours of the unit, (b) present working hours, (c) timer for 30/60/90/120 min - the unit shall stop by the time over.
- Ultraviolet light in the box on top of the unit to sterilize the used tube.

### SPECIFICATION

Brand	: JAY
Capacity(Litre per minute)	: 5 LPM
Capacity	: 5 L
Operation Mode	: Non-Portable
Condition	: New